

Oclube®

Polyethylene Glycol 400 0.4%
& Propylene Glycol 0.3%

Sterile Eye Drops

Description

Oclube® is the combination of Polyethylene Glycol 400 & Propylene Glycol which is clinically proven to reduce both the measurable signs and patient-described symptoms of dry eye. It provides immediate, long lasting moisture and protection. It prolongs tear film break-up time (TFBUT) with a unique biphasic mechanism of action (MOA). As a result, **Oclube®** offers long-lasting symptom relief and an enhanced environment for ocular surface repair. **Oclube®** eye drops contain a unique gelling and lubricating polymer system formulated to adjust to each patient's individual tear film pH. This ensures optimum performance and provides enhanced relief for a patient's specific physiological needs.

Indications

Oclube® is indicated for the treatment of the temporary relief of burning and irritation due to dryness of the eye.

Dosage & administration

The recommended dose of **Oclube®** is one or two drops in affected eye(s) as needed.

Use in pregnancy and lactation

Pregnancy: Due to the negligible systemic exposure and the lack of pharmacological activity, this medication can be used during pregnancy.

Lactation: It is not known whether this medication is excreted in human milk. Because many drugs are excreted in human milk, caution should be exercised when this product is applied to a lactating woman.

Side effects

Generally well tolerated. A very serious allergic reaction to this drug is rare. It should not be used if allergic condition occurs to any ingredients of the product.

Contraindications

It is contraindicated in known hypersensitivity to any of the component of this product.

Warnings & Precautions

Oclube® is only for topical use and not for injection or oral use. Do not use if this product changes color or becomes cloudy. To prevent contaminating the dropper tip and solution, care should be taken not to touch the eyelid or surrounding areas with the dropper tip of the bottle. Keep bottle tightly closed when not in use.

Pharmaceutical precautions

- Store in a cool (15°C - 25°C) & dry place protected from light.
- Do not use after one month once the cap is opened.

Presentation

Oclube® eye drops: Each ml contains Polyethylene Glycol 400 BP 4 mg & Propylene Glycol BP 3 mg.

Package quantities

Oclube® eye drops: Plastic dropper bottle of 10 ml.

® Registered Trade Mark